



Human PRPF40A blocking peptide (CDBP1533)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-HYPA/FBP11 antibody
Antigen Description	PRPF40A (PRP40 pre-mRNA processing factor 40 homolog A (<i>S. cerevisiae</i>)) is a protein-coding gene. Diseases associated with PRPF40A include rett syndrome, and huntington's disease. GO annotations related to this gene include proline-rich region binding. An important paralog of this gene is PRPF40B.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	PRPF40A PRP40 pre-mRNA processing factor 40 homolog A (<i>S. cerevisiae</i>) [<i>Homo sapiens</i> (human)]
Official Symbol	PRPF40A
Synonyms	PRPF40A; PRP40 pre-mRNA processing factor 40 homolog A (<i>S. cerevisiae</i>); HYPA; FBP11; FLAF1; FNBP3; HIP10; Prp40; FBP-11; HIP-10; NY-REN-6; pre-mRNA-processing factor 40 homolog A; NY-REN-6 antigen; formin binding protein 3; formin-binding protein 11; huntingtin

yeast partner A; Fas-ligand associated factor 1; fas ligand-associated factor 1; Huntingtin-interacting protein A; renal carcinoma antigen NY-REN-6; huntingtin-interacting protein 10;

Entrez Gene ID	55660
mRNA Refseq	NM_017892.3
Protein Refseq	NP_060362.3
UniProt ID	O75400
Chromosome Location	2q23.3
Pathway	Spliceosome, organism-specific biosystem; Spliceosome, conserved biosystem; mRNA processing, organism-specific biosystem;
Function	poly(A) RNA binding; proline-rich region binding; protein binding;