



Human IER5 blocking peptide (CDBP1552)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-IER5 antibody
Antigen Description	This gene encodes a protein that is similar to other immediate early response proteins. In the mouse, a similar gene may play an important role in mediating the cellular response to mitogenic signals. Studies in rats found the expression of a similar gene to be increased after waking and sleep deprivation. [provided by RefSeq, Jul 2008]
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	Ier5 immediate early response 5 [Mus musculus]
Official Symbol	Ier5
Synonyms	IER5; immediate early response 5; immediate early response gene 5 protein;
Entrez Gene ID	15939
mRNA Refseq	NM_010500

