



Human IMP3 blocking peptide (CDBP1597)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-IMP3 antibody
Antigen Description	This gene encodes the human homolog of the yeast Imp3 protein. The protein localizes to the nucleoli and interacts with the U3 snoRNP complex. The protein contains an S4 domain. [provided by RefSeq, Jul 2008]
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	IMP3 IMP3, U3 small nucleolar ribonucleoprotein [Homo sapiens (human)]
Official Symbol	IMP3
Synonyms	IMP3; IMP3, U3 small nucleolar ribonucleoprotein; BRMS2; MRPS4; C15orf12; U3 small nucleolar ribonucleoprotein protein IMP3; U3 snoRNP protein IMP3; U3 snoRNP protein 3 homolog; IMP3, U3 small nucleolar ribonucleoprotein, homolog;
Entrez Gene ID	55272

mRNA Refseq	NM_018285.3
Protein Refseq	NP_060755.1
UniProt ID	Q9NV31
Chromosome Location	15q24
Pathway	Ribosome biogenesis in eukaryotes, organism-specific biosystem; Ribosome biogenesis in eukaryotes, conserved biosystem;
Function	poly(A) RNA binding; protein binding; rRNA binding;