



Human IRAK4 blocking peptide (CDBP1618)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-IRAK4 (N Terminus) antibody
Antigen Description	This gene encodes a kinase that activates NF-kappaB in both the Toll-like receptor (TLR) and T-cell receptor (TCR) signaling pathways. The protein is essential for most innate immune responses. Mutations in this gene result in IRAK4 deficiency and recurrent invasive pneumococcal disease. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2011]
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	IRAK4 interleukin-1 receptor-associated kinase 4 [Homo sapiens]
Official Symbol	IRAK4
Synonyms	IRAK4; interleukin-1 receptor-associated kinase 4; NY REN 64; renal carcinoma antigen NY-REN-64; IPD1; REN64; IRAK-4; NY-REN-64;
Entrez Gene ID	51135

mRNA Refseq	NM_001114182
Protein Refseq	NP_001107654
UniProt ID	Q9NWZ3
Chromosome Location	12q12
Pathway	Activated TLR4 signalling, organism-specific biosystem; Apoptosis, organism-specific biosystem; Apoptosis, conserved biosystem; Chagas disease (American trypanosomiasis), organism-specific biosystem; Chagas disease (American trypanosomiasis), conserved biosystem; Cytokine Signaling in Immune system, organism-specific biosystem; IL-1 Signaling Pathway, organism-specific biosystem;
Function	ATP binding; interleukin-1 receptor binding; magnesium ion binding; nucleotide binding; protein binding; protein kinase activity; protein serine/threonine kinase activity; protein serine/threonine kinase activity;