



Human KIF6 blocking peptide (CDBP1686)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-KIF6 antibody
Antigen Description	KIF6 (kinesin family member 6) is a protein-coding gene. Diseases associated with KIF6 include coronary heart disease, and ischemic heart disease. GO annotations related to this gene include microtubule motor activity and ATP binding. An important paralog of this gene is KIF9.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	KIF6 kinesin family member 6 [Homo sapiens]
Official Symbol	KIF6
Synonyms	KIF6; kinesin family member 6; C6orf102, chromosome 6 open reading frame 102; kinesin-like protein KIF6; dJ137F1.4; dJ188D3.1; dJ1043E3.1; DKFZp451I2418; MGC33317; C6orf102;
Entrez Gene ID	221458

mRNA Refseq	NM_145027
Protein Refseq	NP_659464
UniProt ID	Q6ZMV9
Chromosome Location	6p21.2
Function	ATP binding; microtubule motor activity; nucleotide binding; protein binding;