



Human KLHL15 blocking peptide (CDBP1696)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-KLHL15 antibody
Antigen Description	This gene encodes a member of the kelch-like family of proteins that share a common domain structure consisting of an N-terminal broad-complex, tramtrack, bric-a-brac/poxvirus and zinc finger domain and C-terminal kelch repeat motifs. The encoded protein may be involved in protein ubiquitination and cytoskeletal organization. (provided by RefSeq, Apr 2009)
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	KLHL15 kelch-like 15 (Drosophila) [Homo sapiens]
Official Symbol	KLHL15
Synonyms	KLHL15; kelch-like 15 (Drosophila); kelch-like protein 15; KIAA1677; MGC126148;

MGC126149;

Entrez Gene ID	80311
mRNA Refseq	NM_030624
Protein Refseq	NP_085127
UniProt ID	Q96M94
Chromosome Location	Xp22.1-p21