



Human ANKRD54 blocking peptide (CDBP1748)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-LIAR antibody
Antigen Description	ANKRD54 (ankyrin repeat domain 54) is a protein-coding gene. Diseases associated with ANKRD54 include root resorption, and dental caries. GO annotations related to this gene include protein complex binding and protein kinase regulator activity. An important paralog of this gene is PSMD10.
Species	Human
Conjugate	Unconjugated
Applications	BL, WB
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	ANKRD54 ankyrin repeat domain 54 [Homo sapiens (human)]
Official Symbol	ANKRD54

Synonyms	ANKRD54; ankyrin repeat domain 54; LIAR; ankyrin repeat domain-containing protein 54; lyn-interacting ankyrin repeat protein;
Entrez Gene ID	129138
mRNA Refseq	NM_138797.2
Protein Refseq	NP_620152.1
UniProt ID	Q6NXT1
Chromosome Location	22q13.1
Function	protein complex binding; protein kinase regulator activity;