



Human STK11IP blocking peptide (CDBP1756)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-LIP1/LKB1IP antibody
Antigen Description	STK11IP (serine/threonine kinase 11 interacting protein) is a protein-coding gene. Diseases associated with STK11IP include bannayan-riley-ruvalcaba syndrome, and ruvalcaba syndrome. GO annotations related to this gene include protein kinase binding. An important paralog of this gene is NISCH.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	STK11IP serine/threonine kinase 11 interacting protein [Homo sapiens]
Official Symbol	STK11IP
Synonyms	STK11IP; serine/threonine kinase 11 interacting protein; serine/threonine-protein kinase 11-interacting protein; KIAA1898; LIP1; LKB1 interacting protein; LKB1IP; STK11IP1; STK11 interacting protein; LKB1-interacting protein 1; serine/threonine kinase 11-interacting protein;
Entrez Gene ID	114790

mRNA Refseq	NM_052902
Protein Refseq	NP_443134
UniProt ID	Q8N1F8
Chromosome Location	2q35
Pathway	LKB1 signaling events, organism-specific biosystem; TGF-beta Receptor Signaling Pathway, organism-specific biosystem;
Function	protein kinase binding;