



Human LNP1 blocking peptide (CDBP1773)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-LNP1 antibody
Antigen Description	LNP1 (leukemia NUP98 fusion partner 1) is a protein-coding gene. Diseases associated with LNP1 include leukemia, and acute lymphoblastic leukemia.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	LNP1 leukemia NUP98 fusion partner 1 [Homo sapiens]
Official Symbol	LNP1
Synonyms	LNP1; leukemia NUP98 fusion partner 1;
Entrez Gene ID	348801
mRNA Refseq	NM_001085451
Protein Refseq	NP_001078920

UniProt ID	A1A4G5
-------------------	--------

Chromosome Location	3q12.2
----------------------------	--------