



Human LRFN2 blocking peptide (CDBP1777)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-LRFN2 antibody
Antigen Description	LRFN2 (leucine rich repeat and fibronectin type III domain containing 2) is a protein-coding gene. Diseases associated with LRFN2 include esophageal squamous cell carcinoma, and esophagitis. An important paralog of this gene is LRFN3.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	LRFN2 leucine rich repeat and fibronectin type III domain containing 2 [Homo sapiens (human)]
Official Symbol	LRFN2
Synonyms	LRFN2; leucine rich repeat and fibronectin type III domain containing 2; SALM1; FIGLER2;

KIAA1246; leucine-rich repeat and fibronectin type-III domain-containing protein 2; synaptic adhesion-like molecule 1; fibronectin type III, immunoglobulin and leucine rich repeat domains 2;

Entrez Gene ID	57497
----------------	-----------------------

mRNA Refseq	NM_020737.1
-------------	-----------------------------

Protein Refseq	NP_065788.1
----------------	-----------------------------

UniProt ID	Q9ULH4
------------	--------

Chromosome Location	6p21.2-p21.1
---------------------	--------------
