



Human LRFN3 blocking peptide (CDBP1778)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-LRFN3 antibody
Antigen Description	LRFN3 (leucine rich repeat and fibronectin type III domain containing 3) is a protein-coding gene. Diseases associated with LRFN3 include neuronitis. An important paralog of this gene is LRRN3.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	LRFN3 leucine rich repeat and fibronectin type III domain containing 3 [Homo sapiens]
Official Symbol	LRFN3
Synonyms	LRFN3; leucine rich repeat and fibronectin type III domain containing 3; leucine-rich repeat and fibronectin type-III domain-containing protein 3; fibronectin type III; immunoglobulin and leucine

rich repeat domains 1; FIGLER1; MGC2656; SALM4; synaptic adhesion-like molecule 4;
fibronectin type III, immunoglobulin and leucine rich repeat domains 1;

Entrez Gene ID	79414
mRNA Refseq	NM_024509
Protein Refseq	NP_078785
UniProt ID	Q9BTN0
Chromosome Location	19q13.13