



Human LRFN4 blocking peptide (CDBP1779)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-LRFN4 antibody
Antigen Description	LRFN4 (leucine rich repeat and fibronectin type III domain containing 4) is a protein-coding gene. Diseases associated with LRFN4 include bipolar disorder, and neuronitis. An important paralog of this gene is LRFN3.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	LRFN4 leucine rich repeat and fibronectin type III domain containing 4 [Homo sapiens (human)]
Official Symbol	LRFN4
Synonyms	LRFN4; leucine rich repeat and fibronectin type III domain containing 4; SALM3; SALM3.;

FIGLER6; leucine-rich repeat and fibronectin type-III domain-containing protein 4; fibronectin type III, immunoglobulin and leucine rich repeat domains 6;

Entrez Gene ID	78999
mRNA Refseq	NM_024036.4
Protein Refseq	NP_076941.2
UniProt ID	Q6PJG9
Chromosome Location	11q13.2