



Human LRFN5 blocking peptide (CDBP1780)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-LRFN5 antibody
Antigen Description	LRFN5 (leucine rich repeat and fibronectin type III domain containing 5) is a protein-coding gene. Diseases associated with LRFN5 include neuronitis. An important paralog of this gene is LRFN3.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	LRFN5 leucine rich repeat and fibronectin type III domain containing 5 [Homo sapiens]
Official Symbol	LRFN5
Synonyms	LRFN5; leucine rich repeat and fibronectin type III domain containing 5; C14orf146, chromosome 14 open reading frame 146; leucine-rich repeat and fibronectin type-III domain-

containing protein 5; fibronectin type III; immunoglobulin and leucine rich repeat domains 8; FIGLER8; SALM5; fibronectin type III, immunoglobulin and leucine rich repeat domains 8; C14orf146; FLJ30803; DKFZp686G0210;

Entrez Gene ID	145581
mRNA Refseq	NM_152447
Protein Refseq	NP_689660
UniProt ID	Q96NI6
Chromosome Location	14q21.1