



Human LRRC19 blocking peptide (CDBP1784)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-LRRC19 antibody
Antigen Description	LRRC19 (leucine rich repeat containing 19) is a protein-coding gene. An important paralog of this gene is SLITRK2.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	LRRC19 leucine rich repeat containing 19 [Homo sapiens]
Official Symbol	LRRC19
Synonyms	LRRC19; leucine rich repeat containing 19; leucine-rich repeat-containing protein 19; FLJ21302;

Entrez Gene ID	64922
mRNA Refseq	NM_022901
Protein Refseq	NP_075052
UniProt ID	Q9H756
Chromosome Location	9p21.1