



# Mouse LY6E blocking peptide (CDBP1792)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                            |                                                                                                                                                                             |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Overview</b>    | Blocking/Immunizing peptide for anti-LY6E (mouse) antibody                                                                                                                  |
| <b>Antigen Description</b> | LY6E (lymphocyte antigen 6 complex, locus E) is a protein-coding gene. Diseases associated with LY6E include subvalvular aortic stenosis, and acute promyelocytic leukemia. |
| <b>Species</b>             | Mouse                                                                                                                                                                       |
| <b>Conjugate</b>           | Unconjugated                                                                                                                                                                |
| <b>Applications</b>        | Apuri, BL, ELISA                                                                                                                                                            |
| <b>Format</b>              | Lyophilized powder                                                                                                                                                          |
| <b>Size</b>                | 100 µg                                                                                                                                                                      |
| <b>Preservative</b>        | None                                                                                                                                                                        |
| <b>Storage</b>             | Shipped at ambient temperature, store at -20°C.                                                                                                                             |

## GENE INFORMATION

|                        |                                                                                                                                                                 |
|------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>       | <a href="#">Ly6e lymphocyte antigen 6 complex, locus E [ Mus musculus ]</a>                                                                                     |
| <b>Official Symbol</b> | LY6E                                                                                                                                                            |
| <b>Synonyms</b>        | LY6E; lymphocyte antigen 6 complex, locus E; lymphocyte antigen 6E; ly-6E; stem cell antigen 2; thymic shared antigen 1; 9804; Ly67; Tsa1; RIG-E; Sca-2; TSA-1; |
| <b>Entrez Gene ID</b>  | <a href="#">17069</a>                                                                                                                                           |
| <b>mRNA Refseq</b>     | <a href="#">NM_001164036</a>                                                                                                                                    |

