



Human LYPLAL1 blocking peptide (CDBP1794)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-LYPLAL1 antibody
Antigen Description	LYPLAL1 (lysophospholipase-like 1) is a protein-coding gene. Diseases associated with LYPLAL1 include peters anomaly, and obesity. GO annotations related to this gene include lysophospholipase activity. An important paralog of this gene is LYPLA2.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	LYPLAL1 lysophospholipase-like 1 [Homo sapiens]
Official Symbol	LYPLAL1
Synonyms	LYPLAL1; lysophospholipase-like 1; lysophospholipase-like protein 1; Q96AV0; FLJ99730; KIAA1238;
Entrez Gene ID	127018
mRNA Refseq	NM_138794

Protein Refseq	NP_620149
UniProt ID	Q5VWZ2
Chromosome Location	1q41
Function	hydrolase activity; lysophospholipase activity;