



Human LYRM2 blocking peptide (CDBP1795)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-LYRM2 antibody
Antigen Description	LYRM2 (LYR motif containing 2) is a protein-coding gene.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	LYRM2 LYR motif containing 2 [Homo sapiens (human)]
Official Symbol	LYRM2
Synonyms	LYRM2; LYR motif containing 2; DJ122O8.2; LYR motif-containing protein 2;
Entrez Gene ID	57226

mRNA Refseq	NM_020466.4
Protein Refseq	NP_065199.1
UniProt ID	Q9NU23
Chromosome Location	6q15