



Human NAA35 blocking peptide (CDBP1809)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-MAK10 antibody
Antigen Description	NAA35 (N(alpha)-acetyltransferase 35, NatC auxiliary subunit) is a protein-coding gene. Diseases associated with NAA35 include malaria.
Species	Human
Conjugate	Unconjugated
Applications	BL, WB
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	NAA35 N(alpha)-acetyltransferase 35, NatC auxiliary subunit [Homo sapiens (human)]
Official Symbol	NAA35
Synonyms	NAA35; N(alpha)-acetyltransferase 35, NatC auxiliary subunit; EGAP; MAK10; MAK10P; bA379P1.1; N-alpha-acetyltransferase 35, NatC auxiliary subunit; protein MAK10 homolog; corneal wound healing-related protein; embryonic growth-associated protein homolog; MAK10

homolog, amino-acid N-acetyltransferase subunit;

Entrez Gene ID	60560
mRNA Refseq	NM_024635.3
Protein Refseq	NP_078911.3
UniProt ID	Q5VZE5
Chromosome Location	9q21.33
Function	protein binding;