



Human MAPRE3 blocking peptide (CDBP1819)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-MAPRE3 antibody
Antigen Description	The protein encoded by this gene is a member of the RP/EB family of genes. The protein localizes to the cytoplasmic microtubule network and binds APCL, a homolog of the adenomatous polyposis coli tumor suppressor gene. [provided by RefSeq, Jul 2008]
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	MAPRE3 microtubule-associated protein, RP/EB family, member 3 [Homo sapiens (human)]
Official Symbol	MAPRE3
Synonyms	MAPRE3; microtubule-associated protein, RP/EB family, member 3; EB3; RP3; EBF3; EBF3-S; microtubule-associated protein RP/EB family member 3; APC binding protein; end-binding protein 3; EB1 protein family member 3;
Entrez Gene ID	22924

mRNA Refseq	NM_012326.2
Protein Refseq	NP_036458.2
UniProt ID	Q9UPY8
Chromosome Location	2p23.3-p23.1
Function	microtubule binding; protein binding;