



Human LY86 blocking peptide (CDBP1845)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-MD-1 antibody
Antigen Description	LY86 (lymphocyte antigen 86) is a protein-coding gene. Diseases associated with LY86 include kawasaki disease, and tonsillitis.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	LY86 lymphocyte antigen 86 [Homo sapiens (human)]
Official Symbol	LY86
Synonyms	LY86; lymphocyte antigen 86; MD-1; MMD-1; dJ80N2.1; ly-86; protein MD-1; MD-1, RP105-associated;

Entrez Gene ID	9450
mRNA Refseq	NM_004271.3
Protein Refseq	NP_004262.1
UniProt ID	O95711
Chromosome Location	6p25.1
Pathway	Immune System, organism-specific biosystem; Innate Immune System, organism-specific biosystem; Toll Like Receptor 4 (TLR4) Cascade, organism-specific biosystem; Toll-Like Receptors Cascades, organism-specific biosystem;