



Human METTL7A blocking peptide (CDBP1866)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-MettL7A antibody
Antigen Description	METTL7A (methyltransferase like 7A) is a protein-coding gene. Diseases associated with METTL7A include tuberculosis. GO annotations related to this gene include methyltransferase activity. An important paralog of this gene is METTL7B.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	METTL7A methyltransferase like 7A [Homo sapiens]
Official Symbol	METTL7A
Synonyms	METTL7A; methyltransferase like 7A; methyltransferase-like protein 7A; DKFZP586A0522; protein AAM-B; AAM-B; DKFZp586A0522;

Entrez Gene ID	25840
mRNA Refseq	NM_014033
Protein Refseq	NP_054752
UniProt ID	Q9H8H3
Chromosome Location	12q13.12
Function	methyltransferase activity; transferase activity;