



Human MFSD2A blocking peptide (CDBP1871)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-MFSD2A antibody
Antigen Description	MFSD2A (major facilitator superfamily domain containing 2A) is a protein-coding gene. Diseases associated with MFSD2A include hepatoblastoma, and choriocarcinoma. An important paralog of this gene is MFSD2B.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	MFSD2A major facilitator superfamily domain containing 2A [Homo sapiens (human)]
Official Symbol	MFSD2A
Synonyms	MFSD2A; major facilitator superfamily domain containing 2A; MFSD2; major facilitator superfamily domain-containing protein 2A;

Entrez Gene ID	84879
mRNA Refseq	NM_001136493.2
Protein Refseq	NP_001129965.1
UniProt ID	Q71RE4
Chromosome Location	1p34.2