



# Human MFSD2B blocking peptide (CDBP1872)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                            |                                                                                                                                    |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Overview</b>    | Blocking peptide for anti-MFSD2B antibody                                                                                          |
| <b>Antigen Description</b> | MFSD2B (major facilitator superfamily domain containing 2B) is a protein-coding gene. An important paralog of this gene is MFSD2A. |
| <b>Species</b>             | Human                                                                                                                              |
| <b>Conjugate</b>           | Unconjugated                                                                                                                       |
| <b>Applications</b>        | BL                                                                                                                                 |
| <b>Format</b>              | Liquid                                                                                                                             |
| <b>Concentration</b>       | 200 µg/ml                                                                                                                          |
| <b>Size</b>                | 50 µg                                                                                                                              |
| <b>Buffer</b>              | PBS containing 0.02% sodium azide                                                                                                  |
| <b>Preservative</b>        | 0.02% Sodium Azide                                                                                                                 |
| <b>Storage</b>             | Store at -20°C, stable for one year.                                                                                               |

## GENE INFORMATION

|                        |                                                                                                                         |
|------------------------|-------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>       | <a href="#">MFSD2B major facilitator superfamily domain containing 2B [ Homo sapiens ]</a>                              |
| <b>Official Symbol</b> | MFSD2B                                                                                                                  |
| <b>Synonyms</b>        | MFSD2B; major facilitator superfamily domain containing 2B; major facilitator superfamily domain-containing protein 2B; |

|                     |                              |
|---------------------|------------------------------|
| Entrez Gene ID      | <a href="#">388931</a>       |
| mRNA Refseq         | <a href="#">NM_001080473</a> |
| Protein Refseq      | <a href="#">NP_001073942</a> |
| UniProt ID          | A6NFX1                       |
| Chromosome Location | 2p23.3                       |