



Human GHITM blocking peptide (CDBP1874)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-MICS1/GHITM (aa117-127) antibody
Antigen Description	GHITM (growth hormone inducible transmembrane protein) is a protein-coding gene. Diseases associated with GHITM include astigmatism, and typhoid fever. GO annotations related to this gene include protein binding and molecular_function.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	GHITM growth hormone inducible transmembrane protein [Homo sapiens]
Official Symbol	GHITM
Synonyms	GHITM; growth hormone inducible transmembrane protein; growth hormone-inducible transmembrane protein; DERP2; HSPC282; My021; PTD010; TMBIM5; transmembrane BAX inhibitor motif containing 5; dermal papilla-derived protein 2; mitochondrial morphology and cristae structure 1; transmembrane BAX inhibitor motif-containing protein 5; MICS1; FLJ26584; DKFZp566C0746;

Entrez Gene ID	27069
mRNA Refseq	NM_014394
Protein Refseq	NP_055209
UniProt ID	Q9H3K2
Chromosome Location	10q23.1
Function	molecular_function;