



Human ZNF667 blocking peptide (CDBP1879)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-MIPU1 antibody
Antigen Description	ZNF667 (zinc finger protein 667) is a protein-coding gene, and is affiliated with the lncRNA class. Among its related super-pathways are Gene Expression. GO annotations related to this gene include DNA binding and zinc ion binding.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	ZNF667 zinc finger protein 667 [Homo sapiens]
Official Symbol	ZNF667
Synonyms	ZNF667; zinc finger protein 667; FLJ14011; myocardial ischemic preconditioning upregulated 1 ortholog; MIPU1; FLJ45518; DKFZp686O111;

Entrez Gene ID	63934
mRNA Refseq	NM_022103
Protein Refseq	NP_071386
UniProt ID	Q5HYK9
Chromosome Location	19q13.43
Function	DNA binding; metal ion binding; zinc ion binding;