



Human MON1A blocking peptide (CDBP1890)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-MON1A (aa478-489) antibody
Antigen Description	MON1A (MON1 secretory trafficking family member A) is a protein-coding gene. Diseases associated with MON1A include space motion sickness, and syringomyelia. An important paralog of this gene is MON1B.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	MON1A MON1 secretory trafficking family member A [Homo sapiens (human)]
Official Symbol	MON1A
Synonyms	MON1A; MON1 secretory trafficking family member A; SAND1; vacuolar fusion protein MON1 homolog A; MON1 homolog A;
Entrez Gene ID	84315
mRNA Refseq	NM_001142501.1

Protein Refseq [NP_001135973.1](#)

UniProt ID Q86VX9

Chromosome Location 3p21.31
