



Human MTBP blocking peptide (CDBP1918)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	MTBP peptide (human)
Antigen Description	This gene encodes a protein that interacts with the oncoprotein mouse double minute 2. The encoded protein regulates progression through the cell cycle and may be involved in tumor formation. (provided by RefSeq, Aug 2012)
Species	Human
Conjugate	Unconjugated
Applications	BL
Concentration	0.2 mg/ml
Size	50 µg
Buffer	Preservative: 0.02% Sodium Azide; Constituents: 0.1% BSA, PBS. pH 7.2
Preservative	0.02% Sodium Azide

GENE INFORMATION

Gene Name	MTBP Mdm2, transformed 3T3 cell double minute 2, p53 binding protein (mouse) binding protein, 104kDa [Homo sapiens]
Official Symbol	MTBP
Synonyms	MTBP; Mdm2, transformed 3T3 cell double minute 2, p53 binding protein (mouse) binding protein, 104kDa; MDM2 (mouse double minute 2) binding protein, 104kD; mdm2-binding protein; hMTBP; MDM2 (mouse double minute 2)-binding protein, 104kD; MDM2BP;
Entrez Gene ID	27085

mRNA Refseq	NM_022045
Protein Refseq	NP_071328
UniProt ID	Q96DY7
Chromosome Location	8q24.1-q24.2