



Human MXRA5 blocking peptide (CDBP1942)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-MXRA5 antibody
Antigen Description	This gene encodes one of the matrix-remodelling associated proteins. This protein contains 7 leucine-rich repeats and 12 immunoglobulin-like C2-type domains related to perlecan. This gene has a pseudogene on chromosome Y. [provided by RefSeq, Mar 2010]
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	MXRA5 matrix-remodelling associated 5 [Homo sapiens]
Official Symbol	MXRA5
Synonyms	MXRA5; matrix-remodelling associated 5; matrix-modeling-associated protein 5; adlican; DKFZp564I1922; adhesion protein with leucine-rich repeats and Immunoglobulin domains related to perlecan;
Entrez Gene ID	25878

mRNA Refseq	NM_015419
Protein Refseq	NP_056234
UniProt ID	Q9NR99
Chromosome Location	Xp22.33