



Human MYRIP blocking peptide (CDBP1952)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-MYRIP/SLAC2C antibody
Antigen Description	MYRIP (myosin VIIA and Rab interacting protein) is a protein-coding gene. Diseases associated with MYRIP include age related macular degeneration, and macular degeneration. GO annotations related to this gene include protein kinase A binding and actin binding. An important paralog of this gene is MLPH.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	MYRIP myosin VIIA and Rab interacting protein [Homo sapiens]
Official Symbol	MYRIP
Synonyms	MYRIP; myosin VIIA and Rab interacting protein; rab effector MyRIP; DKFZp586F1018; exophilin 8; MyRIP; rab effector MYRIP; SLAC2 C; SLAC2C; Slp homologue lacking C2 domains; synaptotagmin like protein homologue lacking C2 domains c; exophilin-8; slp homolog lacking C2 domains c; myosin-VIIa- and Rab-interacting protein; synaptotagmin-like protein

homologue lacking C2 domains-c; SLAC2-C; FLJ44025; MGC130034; MGC130035;

Entrez Gene ID	25924
mRNA Refseq	NM_015460
Protein Refseq	NP_056275
UniProt ID	Q8NFW9
Chromosome Location	3p21.33
Function	Rab GTPase binding; actin binding; metal ion binding; myosin binding; protein binding; zinc ion binding;