



Human EPB41L4A blocking peptide (CDBP1976)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-NBL4 antibody
Antigen Description	Members of the band 4.1 protein superfamily, including EPB41L4A, are thought to regulate the interaction between the cytoskeleton and plasma membrane (Ishiguro et al., 2000 (PubMed 10874211)).
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	EPB41L4A erythrocyte membrane protein band 4.1 like 4A [Homo sapiens]
Official Symbol	EPB41L4A
Synonyms	EPB41L4A; erythrocyte membrane protein band 4.1 like 4A; band 4.1-like protein 4A; NBL4; erythrocyte protein band 4.1-like 4; EPB41L4; FLJ38738; DKFZp566L203;
Entrez Gene ID	64097

mRNA Refseq	NM_022140
Protein Refseq	NP_071423
UniProt ID	Q9HCS5
Chromosome Location	5q21.3
Function	binding; cytoskeletal protein binding; molecular_function;