



Human NTNG1 blocking peptide (CDBP2001)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-Netrin G1 antibody
Antigen Description	Netrin G1 (NTNG1) belongs to a conserved family of proteins that act as axon guidance cues during vertebrate nervous system development (Nakashiba et al., 2000 [PubMed 10964959]).[supplied by OMIM, Mar 2008]
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	NTNG1 netrin G1 [Homo sapiens (human)]
Official Symbol	NTNG1
Synonyms	NTNG1; netrin G1; Lmnt1; netrin-G1; YLSR571; laminet 1; laminet-1; netrin G1f; axon guidance molecule;
Entrez Gene ID	22854
mRNA Refseq	NM_001113226.1

Protein Refseq	NP_001106697.1
UniProt ID	Q5IEC3
Chromosome Location	1p13.3
Pathway	Axon guidance, organism-specific biosystem; Axon guidance, conserved biosystem; Cell adhesion molecules (CAMs), organism-specific biosystem; Cell adhesion molecules (CAMs), conserved biosystem;
Function	protein binding;