



Human NRTN blocking peptide (CDBP2024)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-Neurturin antibody
Antigen Description	Neurturin is a member of the TGF-beta subfamily, TRN. This gene signals through RET and a GPI-linked coreceptor, and promotes survival of neuronal populations. A neurturin mutation has been described in a family with Hirschsprung Disease. [provided by RefSeq, Jul 2008]
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	NRTN neurturin [Homo sapiens (human)]
Official Symbol	NRTN
Synonyms	NRTN; neurturin; NTN; prepro-neurturin;
Entrez Gene ID	4902
mRNA Refseq	NM_004558.3

Protein Refseq	NP_004549.1
UniProt ID	Q99748
Chromosome Location	19p13.3
Pathway	Axon guidance, organism-specific biosystem; Developmental Biology, organism-specific biosystem; NCAM signaling for neurite out-growth, organism-specific biosystem; NCAM1 interactions, organism-specific biosystem;
Function	growth factor activity; receptor binding;