



Human NIPSNAP3A blocking peptide (CDBP2049)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-NIPSNAP3A antibody
Antigen Description	NIPSNAP3A belongs to a family of proteins with putative roles in vesicular transport (Buechler et al., 2004 [PubMed 15177564]).[supplied by OMIM, Mar 2008]
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	NIPSNAP3A nipsnap homolog 3A (C. elegans) [Homo sapiens]
Official Symbol	NIPSNAP3A
Synonyms	NIPSNAP3A; nipsnap homolog 3A (C. elegans); protein NipSnap homolog 3A;

DKFZp564D177; FLJ13953; HSPC299; MGC14553; protein NipSnap homolog 4; target for Salmonella secreted protein C; TASSC; NIPSNAP4;

Entrez Gene ID	25934
mRNA Refseq	NM_015469
Protein Refseq	NP_056284
UniProt ID	Q9UFN0
Chromosome Location	9q31.3
Function	protein binding;