



Human NTN4 blocking peptide (CDBP2093)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-NTN4 (aa 206-217) antibody
Antigen Description	NTN4 belongs to a family of proteins related to laminins (see LAMA1, MIM 150320) Koch et al. (2000) [PubMed 11038171].[supplied by OMIM, Mar 2008]
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	NTN4 netrin 4 [Homo sapiens]
Official Symbol	NTN4
Synonyms	NTN4; netrin 4; netrin-4; beta netrin; Netrin 4; beta-netrin; hepar-derived netrin-like protein; PRO3091; FLJ23180;
Entrez Gene ID	59277
mRNA Refseq	NM_021229

Protein Refseq	NP_067052
UniProt ID	Q9HB63
Chromosome Location	12q22
Pathway	Axon guidance, organism-specific biosystem; Axon guidance, conserved biosystem; Axon guidance, organism-specific biosystem; Developmental Biology, organism-specific biosystem; Netrin-1 signaling, organism-specific biosystem;
Function	laminin-1 binding;