



Human PCNX blocking peptide (CDBP2244)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-Pecanex homolog antibody
Antigen Description	PCNX (pecanex homolog (Drosophila)) is a protein-coding gene. Diseases associated with PCNX include alcohol dependence, and alzheimer's disease. An important paralog of this gene is PCNXL3.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	PCNX pecanex homolog (Drosophila) [Homo sapiens (human)]
Official Symbol	PCNX
Synonyms	PCNX; pecanex homolog (Drosophila); PCNXL1; pecanex-like protein 1; pecanex-like 1;
Entrez Gene ID	22990
mRNA Refseq	NM_014982.2

Protein Refseq	NP_055797.2
UniProt ID	Q96RV3
Chromosome Location	14q24.2