



Mouse PEG12 blocking peptide (CDBP2245)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-Peg12/Frat3 (mouse) antibody
Species	Mouse
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	Peg12 paternally expressed 12 [Mus musculus (house mouse)]
Official Symbol	PEG12
Synonyms	PEG12; paternally expressed 12; Frat3; frequently rearranged in advanced T-cell lymphomas 3;
Entrez Gene ID	27412
mRNA Refseq	NM_013788.2
Protein Refseq	NP_038816.1
UniProt ID	Q9WVA7

Chromosome Location	7 C; 7 34.38 cM
Pathway	Beta-catenin phosphorylation cascade, organism-specific biosystem; Degradation of beta-catenin by the destruction complex, organism-specific biosystem; Signal Transduction, organism-specific biosystem; Signaling by Wnt, organism-specific biosystem;