



Human PGRMC2 blocking peptide (CDBP2273)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-PGRMC2 antibody
Antigen Description	PGRMC2 (progesterone receptor membrane component 2) is a protein-coding gene. GO annotations related to this gene include steroid binding and heme binding. An important paralog of this gene is PGRMC1.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	PGRMC2 progesterone receptor membrane component 2 [Homo sapiens]
Official Symbol	PGRMC2
Synonyms	PGRMC2; progesterone receptor membrane component 2; membrane-associated progesterone receptor component 2; DG6; PMBP; steroid receptor protein DG6; progesterone membrane binding protein;
Entrez Gene ID	10424

mRNA Refseq	NM_006320
Protein Refseq	NP_006311
UniProt ID	O15173
Chromosome Location	4q26
Function	heme binding; receptor activity; steroid binding; steroid hormone receptor activity;