



Human PINX1 blocking peptide (CDBP2307)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-PINX1 (aa225-237) antibody
Antigen Description	PINX1 (PIN2/TERF1 interacting, telomerase inhibitor 1) is a protein-coding gene. Diseases associated with PINX1 include telangiectasis, and lymphedema. GO annotations related to this gene include telomeric RNA binding and telomerase inhibitor activity.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	PINX1 PIN2/TERF1 interacting, telomerase inhibitor 1 [Homo sapiens]
Official Symbol	PINX1
Synonyms	PINX1; PIN2/TERF1 interacting, telomerase inhibitor 1; PIN2/TERF1-interacting telomerase inhibitor 1; FLJ20565; liver related putative tumor suppressor; LPTL; LPTS; MGC8850; PIN2 interacting protein 1; PinX1; 67-11-3 protein; protein 67-11-3; PIN2-interacting protein 1; TRF1-interacting protein 1; pin2-interacting protein X1; liver-related putative tumor suppressor; hepatocellular carcinoma-related putative tumor suppressor;

Entrez Gene ID	54984
mRNA Refseq	NM_017884
Protein Refseq	NP_060354
UniProt ID	Q96BK5
Chromosome Location	8p23
Pathway	Regulation of Telomerase, organism-specific biosystem;
Function	nucleic acid binding; protein binding; telomerase inhibitor activity; telomeric RNA binding;