



Human TMEM115 blocking peptide (CDBP2318)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-PL6/TMEM115 antibody
Antigen Description	TMEM115 (transmembrane protein 115) is a protein-coding gene. Diseases associated with TMEM115 include renal clear cell carcinoma, and malaria.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	TMEM115 transmembrane protein 115 [Homo sapiens (human)]
Official Symbol	TMEM115
Synonyms	TMEM115; transmembrane protein 115; PL6; PP6; placental protein 6;
Entrez Gene ID	11070
mRNA Refseq	NM_007024.4
Protein Refseq	NP_008955.1

UniProt ID	Q12893
Chromosome Location	3p21.3
Function	molecular_function;