



# Human PLAC4 blocking peptide (CDBP2324)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                            |                                                                                                               |
|----------------------------|---------------------------------------------------------------------------------------------------------------|
| <b>Product Overview</b>    | Blocking peptide for anti-PLAC4 antibody                                                                      |
| <b>Antigen Description</b> | PLAC4 (placenta-specific 4) is a protein-coding gene. Diseases associated with PLAC4 include choriocarcinoma. |
| <b>Species</b>             | Human                                                                                                         |
| <b>Conjugate</b>           | Unconjugated                                                                                                  |
| <b>Applications</b>        | BL                                                                                                            |
| <b>Format</b>              | Liquid                                                                                                        |
| <b>Concentration</b>       | 200 µg/ml                                                                                                     |
| <b>Size</b>                | 50 µg                                                                                                         |
| <b>Buffer</b>              | PBS containing 0.02% sodium azide                                                                             |
| <b>Preservative</b>        | 0.02% Sodium Azide                                                                                            |
| <b>Storage</b>             | Store at -20°C, stable for one year.                                                                          |

## GENE INFORMATION

|                        |                                                                                                                                                    |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene Name</b>       | <a href="#">PLAC4 placenta-specific 4 [ Homo sapiens ]</a>                                                                                         |
| <b>Official Symbol</b> | PLAC4                                                                                                                                              |
| <b>Synonyms</b>        | PLAC4; placenta-specific 4; C21orf115, chromosome 21 open reading frame 115; placenta-specific protein 4; MGC126664; MGC126666; PRED78; C21orf115; |

|                     |                           |
|---------------------|---------------------------|
| Entrez Gene ID      | <a href="#">191585</a>    |
| mRNA Refseq         | <a href="#">NM_182832</a> |
| Protein Refseq      | <a href="#">NP_878252</a> |
| UniProt ID          | Q8WY50                    |
| Chromosome Location | 21                        |