



Human DESI2 blocking peptide (CDBP2336)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-PNAS4 antibody
Antigen Description	DESI2 (desumoylating isopeptidase 2) is a protein-coding gene. Diseases associated with DESI2 include lung adenocarcinoma, and adenocarcinoma. GO annotations related to this gene include peptidase activity.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	DESI2 desumoylating isopeptidase 2 [Homo sapiens (human)]
Official Symbol	DESI2
Synonyms	DESI2; desumoylating isopeptidase 2; DESI; DESI1; DeSI-2; PNAS-4; PPPDE1; CGI-146; FAM152A; C1orf121; desumoylating isopeptidase 1; PPPDE peptidase domain containing 1;

PPPDE peptidase domain-containing protein 1; family with sequence similarity 152, member A;

Entrez Gene ID	51029
mRNA Refseq	NM_016076.3
Protein Refseq	NP_057160.2
UniProt ID	Q9BSY9
Chromosome Location	1q44
Function	peptidase activity;