



Human PSPHP1 blocking peptide (CDBP2424)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-PSPHL/CO9 antibody
Antigen Description	This gene is significantly upregulated in Fanconi's anemia fibroblasts but downregulated or absent in fibroblasts from normal donors. It is also highly expressed in FA B-cells of complementation group A. [provided by RefSeq, Jul 2008]
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	PSPHP1 phosphoserine phosphatase pseudogene 1 [Homo sapiens (human)]
Official Symbol	PSPHP1
Synonyms	PSPHP1; phosphoserine phosphatase pseudogene 1; CO9; PSPHL;
Entrez Gene ID	8781
Chromosome Location	7q11.2

Function	phosphoserine phosphatase activity;
-----------------	-------------------------------------
