



Mouse RAB17 blocking peptide (CDBP2450)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-Rab17 (mouse) antibody
Antigen Description	The Rab subfamily of small GTPases plays an important role in the regulation of membrane trafficking. RAB17 is an epithelial cell-specific GTPase (Lutcke et al., 1993 (PubMed 8486736)).(supplied by OMIM, Oct 2009)
Species	Mouse
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	Rab17 RAB17, member RAS oncogene family [Mus musculus]
Official Symbol	RAB17
Synonyms	RAB17; RAB17, member RAS oncogene family; ras-related protein Rab-17; AW413472;
Entrez Gene ID	19329
mRNA Refseq	NM_001159725

Protein Refseq [NP_001153197](#)

Function GTP binding; nucleotide binding;
