



Human RASIP1 blocking peptide (CDBP2484)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-RASIP1/RAIN antibody
Antigen Description	RASIP1 (Ras interacting protein 1) is a protein-coding gene. Diseases associated with RASIP1 include sotos syndrome, and crohn's disease. An important paralog of this gene is RADIL.
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	RASIP1 Ras interacting protein 1 [Homo sapiens (human)]
Official Symbol	RASIP1
Synonyms	RASIP1; Ras interacting protein 1; RAIN; ras-interacting protein 1;
Entrez Gene ID	54922
mRNA Refseq	NM_017805.2
Protein Refseq	NP_060275.2

UniProt ID	Q5U651
Chromosome Location	19q13.33