



Human RASSF4 blocking peptide (CDBP2487)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-RASSF4 antibody
Antigen Description	The function of this gene has not yet been determined but may involve a role in tumor suppression. Alternative splicing of this gene results in several transcript variants; however, most of the variants have not been fully described. [provided by RefSeq, Jul 2008]
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	RASSF4 Ras association (RalGDS/AF-6) domain family member 4 [Homo sapiens]
Official Symbol	RASSF4
Synonyms	RASSF4; Ras association (RalGDS/AF-6) domain family member 4; ras association domain-containing protein 4; AD037; MGC44914; tumor suppressor RASSF4; Ras association domain family 4; Ras association (RalGDS/AF-6) domain family 4;
Entrez Gene ID	83937

mRNA Refseq	NM_032023
Protein Refseq	NP_114412
UniProt ID	Q9H2L5
Chromosome Location	10q11.1
Function	protein binding;