



Human REEP3 blocking peptide (CDBP2498)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-REEP3 antibody
Antigen Description	REEP3 (receptor accessory protein 3) is a protein-coding gene. Diseases associated with REEP3 include alzheimer's disease, and arthritis. An important paralog of this gene is REEP1.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	REEP3 receptor accessory protein 3 [Homo sapiens]
Official Symbol	REEP3
Synonyms	REEP3; receptor accessory protein 3; C10orf74, chromosome 10 open reading frame 74; receptor expression-enhancing protein 3; receptor expression enhancing protein 3; C10orf74; MGC87479;

Entrez Gene ID	221035
mRNA Refseq	NM_001001330
Protein Refseq	NP_001001330
UniProt ID	Q6NUK4
Chromosome Location	10q21.3