



Human RHBDD3 blocking peptide (CDBP2524)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-RHBDD3 antibody
Antigen Description	RHBDD3 (rhomboid domain containing 3) is a protein-coding gene. Diseases associated with RHBDD3 include pituitary tumors. GO annotations related to this gene include serine-type endopeptidase activity.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	RHBDD3 rhomboid domain containing 3 [Homo sapiens]
Official Symbol	RHBDD3
Synonyms	RHBDD3; rhomboid domain containing 3; C22orf3, chromosome 22 open reading frame 3; rhomboid domain-containing protein 3; PTAG; pituitary tumor apoptosis; C22orf3; HS984G1A;

Entrez Gene ID	25807
mRNA Refseq	NM_012265
Protein Refseq	NP_036397
UniProt ID	Q9Y3P4
Chromosome Location	22q12.2
Function	serine-type endopeptidase activity;