



Human RNF31 blocking peptide (CDBP2554)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-RNF31/ZIBRA (N terminus) antibody
Antigen Description	The protein encoded by this gene contains a RING finger, a motif present in a variety of functionally distinct proteins and known to be involved in protein-DNA and protein-protein interactions. [provided by RefSeq, Jul 2008]
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	RNF31 ring finger protein 31 [Homo sapiens]
Official Symbol	RNF31
Synonyms	RNF31; ring finger protein 31; RING finger protein 31; FLJ10111; FLJ23501; ZIBRA; HOIL-1-interacting protein; zinc in-between-RING-finger ubiquitin-associated domain protein; HOIP; MGC19975;
Entrez Gene ID	55072

mRNA Refseq	NM_017999
Protein Refseq	NP_060469
UniProt ID	Q96EP0
Chromosome Location	14q11.2
Function	ligase activity; metal ion binding; protein binding; ubiquitin binding; ubiquitin-protein ligase activity; zinc ion binding;