



Human RNF34 blocking peptide (CDBP2555)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking/Immunizing peptide for anti-RNF34/RFI (internal) antibody
Antigen Description	The protein encoded by this gene contains a RINF finger, a motif known to be involved in protein-protein and protein-DNA interactions. This protein interacts with DNAJA3/hTid-1, which is a DnaJ protein reported to function as a modulator of apoptosis. Overexpression of this gene in Hela cells was shown to confer the resistance to TNF-alpha induced apoptosis, suggesting an anti-apoptotic function of this protein. This protein can be cleaved by caspase-3 during the induction of apoptosis. This protein also targets p53 and phospho-p53 for degradation. Alternatively splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Feb 2012]
Species	Human
Conjugate	Unconjugated
Applications	Apuri, BL, ELISA
Format	Lyophilized powder
Size	100 µg
Preservative	None
Storage	Shipped at ambient temperature, store at -20°C.

GENE INFORMATION

Gene Name	RNF34 ring finger protein 34, E3 ubiquitin protein ligase [Homo sapiens]
Official Symbol	RNF34

Synonyms	RNF34; ring finger protein 34, E3 ubiquitin protein ligase; ring finger protein 34; E3 ubiquitin-protein ligase RNF34; FLJ21786; RIF; RIFF; caspase regulator CARP1; RING finger protein RIFF; FYVE-RING finger protein MOMO; caspases-8 and -10-associated RING finger protein 1; human RING finger homologous to inhibitor of apoptosis protein; RFI; hRFI; CARP1; CARP-1;
Entrez Gene ID	80196
mRNA Refseq	NM_001256858
Protein Refseq	NP_001243787
UniProt ID	Q969K3
Chromosome Location	12q24.31
Pathway	Adaptive Immune System, organism-specific biosystem; Antigen processing: Ubiquitination & Proteasome degradation, organism-specific biosystem; Class I MHC mediated antigen processing & presentation, organism-specific biosystem; Immune System, organism-specific biosystem;
Function	ligase activity; metal ion binding; zinc ion binding;