



Human RPAIN blocking peptide (CDBP2573)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Blocking peptide for anti-RPA Interacting Protein antibody
Antigen Description	RPAIN (RPA interacting protein) is a protein-coding gene, and is affiliated with the lncRNA class. Diseases associated with RPAIN include patent ductus arteriosus, and protein s deficiency. GO annotations related to this gene include protein complex binding and metal ion binding.
Species	Human
Conjugate	Unconjugated
Applications	BL
Format	Liquid
Concentration	200 µg/ml
Size	50 µg
Buffer	PBS containing 0.02% sodium azide
Preservative	0.02% Sodium Azide
Storage	Store at -20°C, stable for one year.

GENE INFORMATION

Gene Name	RPAIN RPA interacting protein [Homo sapiens]
Official Symbol	RPAIN
Synonyms	RPAIN; RPA interacting protein; RPA-interacting protein; hRIP; MGC4189; RIP; nuclear

transporter; RAP interaction protein; HRIP; FLJ25625; FLJ30490; FLJ42429;

Entrez Gene ID	84268
mRNA Refseq	NM_001033002
Protein Refseq	NP_001028174
UniProt ID	Q86UA6
Chromosome Location	17p13.2
Function	metal ion binding; protein complex binding;